

# Mouse Canx/Calnexin ELISA Kit

Catalog No: #EK19262



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Canx/Calnexin ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Calnexin
Target Name	Calnexin
Accession No.	P35564

## Application Details

Assay Type: SandwichDetect Range:0.312-20ng/mL Sensitivity:0.171ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene\_full\_name:Canx

Note: This product is for in vitro research use only and is not intended for use in humans or animals.