

# Chicken CTSL/Cathepsin L1 ELISA Kit

Catalog No: #EK19264



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Chicken CTSL/Cathepsin L1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Chicken
Specificity	Natural and recombinant Chicken Cathepsin L1
Target Name	Cathepsin L1
Other Names	Cathepsin L,CTSL1
Accession No.	P09648

## Application Details

Assay Type: SandwichDetect Range:0.156-10ng/ml Sensitivity:0.063ng/ml Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene\_full\_name:CTSL

Note: This product is for in vitro research use only and is not intended for use in humans or animals.