

Rat Ctsl/Cathepsin L1 ELISA Kit

Catalog No: #EK19267

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Rat Ctsl/Cathepsin L1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Cathepsin L1
Target Name	Cathepsin L1
Other Names	Cathepsin L,Cyclic protein 2,CP-2,Major excreted protein,MEP,Ctsl1
Accession No.	P07154

Application Details

Assay Type: SandwichDetect Range:0.312-20ng/mL Sensitivity:0.156ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:Ctsl

Note: This product is for in vitro research use only and is not intended for use in humans or animals.