

Bovine APOH/Beta-2-glycoprotein 1 ELISA Kit

Catalog No: #EK19270



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Description

Product Name	Bovine APOH/Beta-2-glycoprotein 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Bovine
Specificity	Natural and recombinant Bovine Beta-2-glycoprotein 1
Target Name	Beta-2-glycoprotein 1
Other Names	Apolipoprotein H,Apo-H,Beta-2-glycoprotein I,B2GPI
Accession No.	P17690

Application Details

Assay Type: SandwichDetect Range:3.12-200ng/ml Sensitivity:1.51ng/ml Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene\_full\_name:APOH

Note: This product is for in vitro research use only and is not intended for use in humans or animals.