

# Mouse Apoh/Beta-2-glycoprotein 1 ELISA Kit

Catalog No: #EK19272

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Apoh/Beta-2-glycoprotein 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Beta-2-glycoprotein 1
Target Name	Beta-2-glycoprotein 1
Other Names	APC inhibitor, Activated protein C-binding protein, Apolipoprotein H, Apo-H, Beta-2-glycoprotein I, B2GPI, B2gp1
Accession No.	Q01339

## Application Details

Assay Type: Sandwich Detect Range: 12.5-800ng/ml Sensitivity: 6.28ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: Apoh

Note: This product is for in vitro research use only and is not intended for use in humans or animals.