

Rat Apoh/Beta-2-glycoprotein 1 ELISA Kit

Catalog No: #EK19273

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Rat Apoh/Beta-2-glycoprotein 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Beta-2-glycoprotein 1
Target Name	Beta-2-glycoprotein 1
Other Names	Apolipoprotein H,Apo-H,Beta-2-glycoprotein I,B2GPI
Accession No.	P26644

Application Details

Assay Type: Sandwich Detect Range:1.56-100ng/ml Sensitivity:0.79ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Apoh

Note: This product is for in vitro research use only and is not intended for use in humans or animals.