

Human CDC42/Cell division control protein 42 homolog ELISA Kit

Catalog No: #EK19274

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human CDC42/Cell division control protein 42 homolog ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Cell division control protein 42 homolog
Target Name	Cell division control protein 42 homolog
Other Names	G25K GTP-binding protein
Accession No.	P60953

Application Details

Assay Type: Sandwich Detect Range: 15.6-1000pg/mL Sensitivity: 10pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: CDC42

Note: This product is for in vitro research use only and is not intended for use in humans or animals.