

Rat Cdc42/Cell division control protein 42 homolog
ELISA Kit

Catalog No: #EK19276



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Rat Cdc42/Cell division control protein 42 homolog ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Cell division control protein 42 homolog
Target Name	Cell division control protein 42 homolog
Accession No.	Q8CFN2

Application Details

Assay Type: SandwichDetect Range:15.6-1000pg/mL Sensitivity:7.92pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:Cdc42

Note: This product is for in vitro research use only and is not intended for use in humans or animals.