

Mouse Anxa3/Annexin A3 ELISA Kit

Catalog No: #EK19277



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Anxa3/Annexin A3 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Annexin A3
Target Name	Annexin A3
Other Names	35-alpha calcimedlin,Annexin III,Annexin-3,Lipocortin III,Placental anticoagulant protein III,PAP-III,Anx3
Accession No.	O35639

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.078ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Anxa3

Note: This product is for in vitro research use only and is not intended for use in humans or animals.