

Human C1QC/Complement C1q subcomponent subunit C ELISA Kit

Catalog No: #EK19280



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human C1QC/Complement C1q subcomponent subunit C ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Complement C1q subcomponent subunit C
Target Name	Complement C1q subcomponent subunit C
Other Names	C1QG
Accession No.	P02747

Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/mL Sensitivity: 33pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids. Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: C1QC

Note: This product is for in vitro research use only and is not intended for use in humans or animals.