

# Mouse C1qc/Complement C1q subcomponent subunit C ELISA Kit

Catalog No: #EK19281

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse C1qc/Complement C1q subcomponent subunit C ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Complement C1q subcomponent subunit C
Target Name	Complement C1q subcomponent subunit C
Other Names	C1qg
Accession No.	Q02105

## Application Details

Assay Type: SandwichDetect Range:0.156-10ng/mL Sensitivity:0.082ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene\_full\_name:C1qc

Note: This product is for in vitro research use only and is not intended for use in humans or animals.