

Rat C1qc/Complement C1q subcomponent subunit C ELISA Kit

Catalog No: #EK19282

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Rat C1qc/Complement C1q subcomponent subunit C ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Complement C1q subcomponent subunit C
Target Name	Complement C1q subcomponent subunit C
Other Names	C1qg
Accession No.	P31722

Application Details

Assay Type: SandwichDetect Range:12.5-800ng/mL Sensitivity:10ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:C1qc

Note: This product is for in vitro research use only and is not intended for use in humans or animals.