

Mouse Cpn2/Carboxypeptidase N subunit 2 ELISA Kit

Catalog No: #EK19284



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Mouse Cpn2/Carboxypeptidase N subunit 2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Carboxypeptidase N subunit 2
Target Name	Carboxypeptidase N subunit 2
Other Names	Carboxypeptidase N 83 kDa chain,Carboxypeptidase N large subunit,Carboxypeptidase N polypeptide 2,Carboxypeptidase N regulatory subunit
Accession No.	Q9DBB9

Application Details

Assay Type: SandwichDetect Range:1.56-100ng/mL Sensitivity:0.69ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:Cpn2

Note: This product is for in vitro research use only and is not intended for use in humans or animals.