

# Bovine AVP/Vasopressin-neurophysin 2-copeptin ELISA Kit

Catalog No: #EK19285

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |                                                                   |
|--------------------|-------------------------------------------------------------------|
| Product Name       | Bovine AVP/Vasopressin-neurophysin 2-copeptin ELISA Kit           |
| Brief Description  | ELISA Kits                                                        |
| Applications       | ELISA                                                             |
| Species Reactivity | Bovine                                                            |
| Specificity        | Natural and recombinant Bovine Vasopressin-neurophysin 2-copeptin |
| Target Name        | Vasopressin-neurophysin 2-copeptin                                |
| Other Names        | AVP-NP11                                                          |
| Accession No.      | P01180                                                            |

## Application Details

Assay Type: Sandwich Detect Range: 78-5000 pg/ml Sensitivity: 21 pg/ml Sample Type: Serum, plasma, tissue homogenates and other biological fluids. Sample Volume: 100 µl Assay Time: 3-5 h Detection wavelength: 450 nm gene\_full\_name: AVP

Note: This product is for in vitro research use only and is not intended for use in humans or animals.