

# Mouse Avp/Vasopressin-neurophysin 2-copeptin ELISA Kit

Catalog No: #EK19286



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Avp/Vasopressin-neurophysin 2-copeptin ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Vasopressin-neurophysin 2-copeptin
Target Name	Vasopressin-neurophysin 2-copeptin
Other Names	AVP-NP11
Accession No.	P35455

## Application Details

Assay Type: Sandwich Detect Range: 31.2-2000pg/mL Sensitivity: 15.77pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: Avp

Note: This product is for in vitro research use only and is not intended for use in humans or animals.