

Pig AVP/Vasopressin-neurophysin 2-copeptin ELISA Kit

Catalog No: #EK19287

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Pig AVP/Vasopressin-neurophysin 2-copeptin ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Pig
Specificity	Natural and recombinant Pig Vasopressin-neurophysin 2-copeptin
Target Name	Vasopressin-neurophysin 2-copeptin
Other Names	AVP-NP11
Accession No.	P01183

Application Details

Assay Type: SandwichDetect Range:7.8-500pg/mL Sensitivity:3.7 pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:AVP

Note: This product is for in vitro research use only and is not intended for use in humans or animals.