

Rat Cpa2/Carboxypeptidase A2 ELISA Kit

Catalog No: #EK19290



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Rat Cpa2/Carboxypeptidase A2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Carboxypeptidase A2
Target Name	Carboxypeptidase A2
Accession No.	P19222

Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/mL Sensitivity: 39pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids. Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Cpa2

Note: This product is for in vitro research use only and is not intended for use in humans or animals.