

Rat Cpt1a/Carnitine O-palmitoyltransferase 1, liver isoform ELISA Kit

Catalog No: #EK19292

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Rat Cpt1a/Carnitine O-palmitoyltransferase 1, liver isoform ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Carnitine O-palmitoyltransferase 1, liver isoform
Target Name	Carnitine O-palmitoyltransferase 1, liver isoform
Other Names	CPT1-L, Carnitine O-palmitoyltransferase I, liver isoform, CPT I, Carnitine palmitoyltransferase 1A, Cpt-1, Cpt1
Accession No.	P32198

Application Details

Assay Type: Sandwich Detect Range: 15.6-1000pg/mL Sensitivity: 7.8pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Cpt1a

Note: This product is for in vitro research use only and is not intended for use in humans or animals.