

Rat Actn1/Alpha-actinin-1 ELISA Kit

Catalog No: #EK19294



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Rat Actn1/Alpha-actinin-1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Alpha-actinin-1
Target Name	Alpha-actinin-1
Other Names	Alpha-actinin cytoskeletal isoform,F-actin cross-linking protein,Non-muscle alpha-actinin-1
Accession No.	Q9Z1P2

Application Details

Assay Type: SandwichDetect Range:0.156-10ng/mL Sensitivity:0.078ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:Actn1

Note: This product is for in vitro research use only and is not intended for use in humans or animals.