

# Pig C-Peptide/C-Peptide ELISA Kit

Catalog No: #EK19300



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Pig C-Peptide/C-Peptide ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Pig
Specificity	Natural and recombinant Pig C-Peptide
Target Name	C-Peptide

## Application Details

Assay Type: Sandwich Detect Range: 1.25-80ng/mL Sensitivity: 0.625ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: C-Peptide

Note: This product is for in vitro research use only and is not intended for use in humans or animals.