

Mouse Cd55/Complement decay-accelerating factor, GPI-anchored ELISA Kit

Catalog No: #EK19301



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Cd55/Complement decay-accelerating factor, GPI-anchored ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Complement decay-accelerating factor, GPI-anchored
Target Name	Complement decay-accelerating factor, GPI-anchored
Other Names	DAF-GPI,CD55,Cd55a,Daf,Daf1
Accession No.	Q61475

Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.088ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name: Cd55

Note: This product is for in vitro research use only and is not intended for use in humans or animals.