

Dog BGLAP/Osteocalcin ELISA Kit

Catalog No: #EK19303



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Dog BGLAP/Osteocalcin ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Dog
Specificity	Natural and recombinant Dog Osteocalcin
Target Name	Osteocalcin
Other Names	Bone Gla protein,BGP,Gamma-carboxyglutamic acid-containing protein
Accession No.	P81455

Application Details

Assay Type: Sandwich Detect Range:0.625-40ng/mL Sensitivity:0.317ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:BGLAP

Note: This product is for in vitro research use only and is not intended for use in humans or animals.