

Bovine CSTA/Cystatin-A ELISA Kit

Catalog No: #EK19304



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Bovine CSTA/Cystatin-A ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Bovine
Specificity	Natural and recombinant Bovine Cystatin-A
Target Name	Cystatin-A
Other Names	Cystatin-AS, Stefin-A, STF1
Accession No.	P80416

Application Details

Assay Type: Sandwich Detect Range: 1.56-100ng/mL Sensitivity: 0.82ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: CSTA

Note: This product is for in vitro research use only and is not intended for use in humans or animals.