

# Bovine CKM/Creatine kinase M-type ELISA Kit

Catalog No: #EK19307



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Bovine CKM/Creatine kinase M-type ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Bovine
Specificity	Natural and recombinant Bovine Creatine kinase M-type
Target Name	Creatine kinase M-type
Other Names	Creatine kinase M chain,M-CK
Accession No.	Q9XSC6

## Application Details

Assay Type: Sandwich Detect Range:0.625-40ng/mL Sensitivity:0.38ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:CKM

Note: This product is for in vitro research use only and is not intended for use in humans or animals.