

Dog CKM/Creatine kinase M-type ELISA Kit

Catalog No: #EK19308



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Dog CKM/Creatine kinase M-type ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Dog
Specificity	Natural and recombinant Dog Creatine kinase M-type
Target Name	Creatine kinase M-type
Other Names	Creatine kinase M chain,M-CK
Accession No.	P05123

Application Details

Assay Type: Sandwich Detect Range:0.78-50ng/mL Sensitivity:0.431ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:CKM

Note: This product is for in vitro research use only and is not intended for use in humans or animals.