

Human AMD1/S-adenosylmethionine decarboxylase proenzyme ELISA Kit

Catalog No: #EK19309



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human AMD1/S-adenosylmethionine decarboxylase proenzyme ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human S-adenosylmethionine decarboxylase proenzyme
Target Name	S-adenosylmethionine decarboxylase proenzyme
Other Names	AdoMetDC,AMD
Accession No.	P17707

Application Details

Assay Type: Sandwich Detect Range: 1.56-100ng/mL Sensitivity: 0.83ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: AMD1

Note: This product is for in vitro research use only and is not intended for use in humans or animals.