

# Mouse Amd1/S-adenosylmethionine decarboxylase proenzyme 1 ELISA Kit

Catalog No: #EK19310



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Amd1/S-adenosylmethionine decarboxylase proenzyme 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse S-adenosylmethionine decarboxylase proenzyme 1
Target Name	S-adenosylmethionine decarboxylase proenzyme 1
Other Names	AdoMetDC 1
Accession No.	P0DMN7

## Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/mL Sensitivity: 39.9pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: Amd1

Note: This product is for in vitro research use only and is not intended for use in humans or animals.