

Rat Amd1/S-adenosylmethionine decarboxylase proenzyme ELISA Kit

Catalog No: #EK19311



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Rat Amd1/S-adenosylmethionine decarboxylase proenzyme ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat S-adenosylmethionine decarboxylase proenzyme
Target Name	S-adenosylmethionine decarboxylase proenzyme
Other Names	AdoMetDC
Accession No.	P17708

Application Details

Assay Type: Sandwich Detect Range: 0.312-20ng/mL Sensitivity: 0.156ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Amd1

Note: This product is for in vitro research use only and is not intended for use in humans or animals.