

Mouse Cartpt/Cocaine- and amphetamine-regulated transcript protein ELISA Kit

Catalog No: #EK19313



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Mouse Cartpt/Cocaine- and amphetamine-regulated transcript protein ELISA Kit |
| Brief Description | ELISA Kits |
| Applications | ELISA |
| Species Reactivity | Ms |
| Specificity | Natural and recombinant Mouse Cocaine- and amphetamine-regulated transcript protein |
| Target Name | Cocaine- and amphetamine-regulated transcript protein |
| Other Names | Cart |
| Accession No. | P56388 |

Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.089ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Cartpt

Note: This product is for in vitro research use only and is not intended for use in humans or animals.