

# Bovine AQP4/Aquaporin-4 ELISA Kit

Catalog No: #EK19321



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Bovine AQP4/Aquaporin-4 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Bovine
Specificity	Natural and recombinant Bovine Aquaporin-4
Target Name	Aquaporin-4
Other Names	AQP-4,Mercurial-insensitive water channel,MIWC,WCH4
Accession No.	O77750

## Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.161ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:AQP4

Note: This product is for in vitro research use only and is not intended for use in humans or animals.