

# Pig AQP5/Aquaporin-5 ELISA Kit

Catalog No: #EK19322



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Pig AQP5/Aquaporin-5 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Pig
Specificity	Natural and recombinant Pig Aquaporin-5
Target Name	Aquaporin-5
Other Names	AQP-5
Accession No.	A8W649

## Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.052ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: AQP5

Note: This product is for in vitro research use only and is not intended for use in humans or animals.