

Bovine CCND1/G1/S-specific cyclin-D1 ELISA Kit

Catalog No: #EK19324



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Bovine CCND1/G1/S-specific cyclin-D1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Bovine
Specificity	Natural and recombinant Bovine G1/S-specific cyclin-D1
Target Name	G1/S-specific cyclin-D1
Accession No.	Q2KI22

Application Details

Assay Type: Sandwich Detect Range: 0.312-20ng/mL Sensitivity: 0.165ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: CCND1

Note: This product is for in vitro research use only and is not intended for use in humans or animals.