

# Rat Ccnd2/G1/S-specific cyclin-D2 ELISA Kit

Catalog No: #EK19326



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Rat Ccnd2/G1/S-specific cyclin-D2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat G1/S-specific cyclin-D2
Target Name	G1/S-specific cyclin-D2
Other Names	Vin-1 proto-oncogene, Vin-1
Accession No.	Q04827

## Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.089ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: Ccnd2

Note: This product is for in vitro research use only and is not intended for use in humans or animals.