

# Rat Ccnd3/G1/S-specific cyclin-D3 ELISA Kit

Catalog No: #EK19327



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Rat Ccnd3/G1/S-specific cyclin-D3 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat G1/S-specific cyclin-D3
Target Name	G1/S-specific cyclin-D3
Accession No.	P48961

## Application Details

Assay Type: Sandwich Detect Range: 3.12-200ng/mL Sensitivity: 1.63ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: Ccnd3

Note: This product is for in vitro research use only and is not intended for use in humans or animals.