

Rat Cpb2/Carboxypeptidase B2 ELISA Kit

Catalog No: #EK19329



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Rat Cpb2/Carboxypeptidase B2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Carboxypeptidase B2
Target Name	Carboxypeptidase B2
Other Names	Carboxypeptidase R,CPR,Carboxypeptidase U,CPU,Thrombin-activable fibrinolysis inhibitor,TAFI
Accession No.	Q9EQV9

Application Details

Assay Type: Sandwich Detect Range:78-5000pg/mL Sensitivity:39pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids. Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Cpb2

Note: This product is for in vitro research use only and is not intended for use in humans or animals.