

Human ALDOB/Fructose-bisphosphate aldolase B ELISA Kit

Catalog No: #EK19330



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|----------------------------------------------------------------|
| Product Name | Human ALDOB/Fructose-bisphosphate aldolase B ELISA Kit |
| Brief Description | ELISA Kits |
| Applications | ELISA |
| Species Reactivity | Hu |
| Specificity | Natural and recombinant Human Fructose-bisphosphate aldolase B |
| Target Name | Fructose-bisphosphate aldolase B |
| Other Names | Liver-type aldolase,ALDB |
| Accession No. | P05062 |

Application Details

Assay Type: Sandwich Detect Range:78-5000pg/mL Sensitivity:39.5pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:ALDOB

Note: This product is for in vitro research use only and is not intended for use in humans or animals.