

Mouse Aldob/Fructose-bisphosphate aldolase B ELISA Kit

Catalog No: #EK19331



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Aldob/Fructose-bisphosphate aldolase B ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Fructose-bisphosphate aldolase B
Target Name	Fructose-bisphosphate aldolase B
Other Names	Aldolase 2,Liver-type aldolase,Aldo2
Accession No.	Q91Y97

Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.12ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Aldob

Note: This product is for in vitro research use only and is not intended for use in humans or animals.