

# Rat Aldob/Fructose-bisphosphate aldolase B ELISA Kit

Catalog No: #EK19332



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Rat Aldob/Fructose-bisphosphate aldolase B ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Fructose-bisphosphate aldolase B
Target Name	Fructose-bisphosphate aldolase B
Other Names	Liver-type aldolase
Accession No.	P00884

## Application Details

Assay Type: Sandwich Detect Range: 0.78-50ng/mL Sensitivity: 0.443ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: Aldob

Note: This product is for in vitro research use only and is not intended for use in humans or animals.