

Mouse Bid/BH3-interacting domain death agonist ELISA Kit

Catalog No: #EK19334



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Bid/BH3-interacting domain death agonist ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse BH3-interacting domain death agonist
Target Name	BH3-interacting domain death agonist
Other Names	p22 BID,BID
Accession No.	P70444

Application Details

Assay Type: Sandwich Detect Range: 1.56-100ng/mL Sensitivity: 0.78ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Bid

Note: This product is for in vitro research use only and is not intended for use in humans or animals.