

# Bovine ANGPT4/Angiopoietin-4 ELISA Kit

Catalog No: #EK19346



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Bovine ANGPT4/Angiopoietin-4 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Bovine
Specificity	Natural and recombinant Bovine Angiopoietin-4
Target Name	Angiopoietin-4
Other Names	ANG-4
Accession No.	Q24K15

## Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.157ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:ANGPT4

Note: This product is for in vitro research use only and is not intended for use in humans or animals.