

# Mouse Cd36/Platelet glycoprotein 4 ELISA Kit

Catalog No: #EK19348



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Cd36/Platelet glycoprotein 4 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Platelet glycoprotein 4
Target Name	Platelet glycoprotein 4
Other Names	Glycoprotein IIIb,GPIIIB,PAS IV,PAS-4,Platelet glycoprotein IV,GPIV,CD36
Accession No.	Q08857

## Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.082ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name: Cd36

Note: This product is for in vitro research use only and is not intended for use in humans or animals.