

Rat Cd36/Platelet glycoprotein 4 ELISA Kit

Catalog No: #EK19349

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Cd36/Platelet glycoprotein 4 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Platelet glycoprotein 4
Target Name	Platelet glycoprotein 4
Other Names	Adipocyte membrane protein,Fatty acid translocase,Fatty acid transport protein,Glycoprotein IIIb,GPIIIB,PAS IV,PAS-4,Platelet glycoprotein IV,GPIV,CD36,Fat
Accession No.	Q07969

Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.156ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name: Cd36

Note: This product is for in vitro research use only and is not intended for use in humans or animals.