

Rat Aldh1a1/Retinal dehydrogenase 1 ELISA Kit

Catalog No: #EK19350

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Aldh1a1/Retinal dehydrogenase 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Retinal dehydrogenase 1
Target Name	Retinal dehydrogenase 1
Other Names	RALDH 1,ALDH-E1,ALHDII,Aldehyde dehydrogenase family 1 member A1,Aldehyde dehydrogenase, cytosolic,Aldh
Accession No.	P51647

Application Details

Assay Type: Sandwich Detect Range:0.625-40ng/ml Sensitivity:0.328ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Aldh1a1

Note: This product is for in vitro research use only and is not intended for use in humans or animals.