

Mouse Adh1/Alcohol dehydrogenase 1 ELISA Kit

Catalog No: #EK19351



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Adh1/Alcohol dehydrogenase 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Alcohol dehydrogenase 1
Target Name	Alcohol dehydrogenase 1
Other Names	ADH-A2,Alcohol dehydrogenase A subunit,Adh-1
Accession No.	P00329

Application Details

Assay Type: Sandwich Detect Range:78-5000pg/mL Sensitivity:39.1pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Adh1

Note: This product is for in vitro research use only and is not intended for use in humans or animals.