

Rat Calb2/Calretinin ELISA Kit

Catalog No: #EK19353



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Rat Calb2/Calretinin ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Calretinin
Target Name	Calretinin
Other Names	CR
Accession No.	P47728

Application Details

Assay Type: Sandwich Detect Range: 15.6-1000pg/mL Sensitivity: 7.98pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Calb2

Note: This product is for in vitro research use only and is not intended for use in humans or animals.