

# Mouse Ccl24/C-C motif chemokine 24 ELISA Kit

Catalog No: #EK19355



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Ccl24/C-C motif chemokine 24 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse C-C motif chemokine 24
Target Name	C-C motif chemokine 24
Other Names	Eosinophil chemotactic protein 2,Eotaxin-2,Small-inducible cytokine A24,Scya24
Accession No.	Q9JKC0

## Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/ml Sensitivity:0.085ng/ml Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:Ccl24

Note: This product is for in vitro research use only and is not intended for use in humans or animals.