

Bovine CEBPB/CCAAT/enhancer-binding protein beta ELISA Kit

Catalog No: #EK19359



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Bovine CEBPB/CCAAT/enhancer-binding protein beta ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Bovine
Specificity	Natural and recombinant Bovine CCAAT/enhancer-binding protein beta
Target Name	CCAAT/enhancer-binding protein beta
Other Names	C/EBP beta
Accession No.	O02755

Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.064ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: CEBPB

Note: This product is for in vitro research use only and is not intended for use in humans or animals.