

Human CEBPB/CCAAT/enhancer-binding protein beta ELISA Kit

Catalog No: #EK19360



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human CEBPB/CCAAT/enhancer-binding protein beta ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human CCAAT/enhancer-binding protein beta
Target Name	CCAAT/enhancer-binding protein beta
Other Names	C/EBP beta,Liver activator protein,LAP,Liver-enriched inhibitory protein,LIP,Nuclear factor NF-IL6,Transcription factor 5,TCF-5,TCF5,PP9092
Accession No.	P17676

Application Details

Assay Type: SandwichDetect Range:0.78-50ng/mL Sensitivity:0.39ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:CEBPB

Note: This product is for in vitro research use only and is not intended for use in humans or animals.