

Rat Bace1/Beta-secretase 1 ELISA Kit

Catalog No: #EK19362



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Rat Bace1/Beta-secretase 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Beta-secretase 1
Target Name	Beta-secretase 1
Other Names	Aspartyl protease 2,ASP2,Beta-site amyloid precursor protein cleaving enzyme 1,Beta-site APP cleaving enzyme 1,Memapsin-2,Membrane-associated aspartic protease 2,Bace
Accession No.	P56819

Application Details

Assay Type: Sandwich Detect Range: 1.56-100ng/mL Sensitivity: 0.71ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Bace1

Note: This product is for in vitro research use only and is not intended for use in humans or animals.