

Mouse Cebpe/CCAAT/enhancer-binding protein epsilon ELISA Kit

Catalog No: #EK19363



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Cebpe/CCAAT/enhancer-binding protein epsilon ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse CCAAT/enhancer-binding protein epsilon
Target Name	CCAAT/enhancer-binding protein epsilon
Other Names	C/EBP epsilon,Gm294
Accession No.	Q6PZD9

Application Details

Assay Type: Sandwich Detect Range:0.78-50ng/mL Sensitivity:0.195ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Cebpe

Note: This product is for in vitro research use only and is not intended for use in humans or animals.